

DATA SHEET TASKFORCE – First Aid & Bleed Control Trauma Kit

SKU Number **SS-6010125060**

Date of Publication **01-09-2025**

EAN13 Barcode **5060660850907**

Commodity Code **30065000.00**

Product Image



Product Description

TASKFORCE – First Aid & Bleed Control Trauma Kit (complies with BS 8599-1:2019)

When severe haemorrhage control is needed, having a kit equipped with the right medical products is essential. Simple, user-friendly Bleed Control Kits are critical to the effectiveness of first responders.

The **SOVOS TASKFORCE** is a hybrid kit for First Aid and Bleed Control Trauma Kit. Combining the elements of a Bleed Control Kitt with that of a Small Workplace First Aid Kit complying with BS 8599-1:2019; this kit is packed with high-quality medical supplies necessary to manage life-threatening haemorrhages, lacerations, burns, and minor injuries.

This kit allows trained First Aiders, Advanced Professionals, and First Responders to provide lifesaving care in critical situations whilst waiting for Emergency Services to arrive.

Contained in a durable weather resistant bag made from advanced materials, this kit is designed for quick access in emergencies. The elements of the First Aid and Bleed Control are separated into two separate internal bags distinguished by colour.

Contents	<u>Main Bag:</u>		
	1	x	Tactical Multi Bandage
	1	x	SAM SPLINT 9" Flat
	1	x	High Performance Emergency Thermal Blanket
	3	x	Security Seals
	<u>Green Bag (Fist Aid):</u>		
	Small Workplace First Aid Kit conforming to BS 8599-1:2019		
	<u>Red Bag (Bleed Control):</u>		
	1	x	SAM XT Tourniquet
	1	x	SWAT-T Rescue Orange Tourniquet
	1	x	6" FCP Military Emergency Bandage
	1	x	ChitoSAM Haemostatic Gauze 4' Z-Fold
	1	x	Sterile Saline Spray (100ml)
	1	x	High Performance Medical Shears
	1	x	Medical Marker Pen
2	x	Disposable Medical Gloves	
2	x	De-Contamination Wipes 30cm x 20cm	

Safety Standards	All medical products used in the SOVOS Kits are CE Marked and conform with the relevant European Union harmonisation legislation.
------------------	---

Product Weight	2300g
Size	320mm x 270mm x 120mm (approx.)

Storage	- Keep stored and transported at ambient temperature (0°C to + 30°C) - Do not store in direct/high heat or sunlight as this may affect the internal components of the kit
Lifetime	- Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. - Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable.
Cleaning	Clean the storage bag using warm soapy water and soft cloth
Disposal	Recycling is preferred where possible. Look out for the recycling symbol for material categories

Replacement Components	• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents • Refer to the SOVOS Website for replacement items www.sovoshelmets.com
Additional Information	• In extreme cases of haemorrhage further products may be required over and above the items in these kits. • This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019
Training	• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use. • Training Providers can be recommended upon request
Trademarks	• SOVOS® is a registered trademark